

# DarkAdapt

Dark adaptation made easy

Striatech's dark adaptation box takes the hassle out of dark-adapting your experimental animals. DarkAdapt is completely light-tight. Nevertheless, it provides robust ventilation, so that you can thoroughly dark adapt your animal even over extended periods of time, and even in a well-lit laboratory environment.



**STRIA.TECH**

Your Experts for Neuroscience Testing Tools

[www.stria.tech](http://www.stria.tech)



# STRIA.TECH

Your Experts for Neuroscience Testing Tools

## Technical specifications

- ▶ **Absolutely light tight**
- ▶ **Robust ventilation**
  - Ultra low-noise fan
  - More than 100-fold air exchange per hour
- ▶ **Size**
  - Inside: 30 cm x 53 cm x 29 cm
  - Outside: 35 cm x 55 cm x 35 cm
  - Easily fits Type II-long mouse cage or rat cage
- ▶ **Stationary and mobile use**
  - Keep in lab
  - Transport between rooms or facilities

**Contact**    **Striatech GmbH**  
Vor dem Kreuzberg 17  
72070 Tübingen  
Germany

 +49 7071 53913-0

 [info@stria.tech](mailto:info@stria.tech)

 [www.stria.tech](http://www.stria.tech)

Developed with support from:

Gefördert durch:

